

APPENDIX 5. LIST OF POST-INTERVENTION VERIFICATION

(mark with an "X" while the steps are being completed in the following order)	
1. Remove and dispose of PPE in the same room where the procedure is performed.	<input type="checkbox"/>
<p>It is useful to follow a visual guide during the removal of the different elements, such as the one provided by the WHO: https://www.who.int/csr/resources/publications/ebola/ppe-steps/en/</p>	
2. Dispose of the used one-use materials into a container for infectious biological waste (Category B, UN 3291).	<input type="checkbox"/>
3. Disconnect the pH/Impedance tube from the recorder before removing it to facilitate its fast disposal into a plastic bag.	<input type="checkbox"/>
4. Check in perfusion systems: do not close the perfusion until the catheter has been removed from the patient.	<input type="checkbox"/>
5. Disinfection with a viricidal agent of stretcher, railing, and other elements with which the patient or the physician has had contact.	<input type="checkbox"/>
6. pH or pH-impedance recorder, cover, and security strap: place them in a container and disinfect with 70% alcohol or other viricidal agent before performing the data transference.	<input type="checkbox"/>
7. Disinfection of the rest of the devices according to the usual protocol of the Motility Unit.	<input type="checkbox"/>